**VALIVACQ software**

Available languages : French, English, Italian, German, Spanish

**ValiVACQ has been designed for validation of hospital sterilization. EN554 compatible.**


ValiVACQ is a validation software enabling the edition of reports in compliance with EN554 norm.

User friendly and simple of use, it allows a precise analysis of the physical magnitudes measured during the whole cycle (**temps d'équilibrage, temps de maintien...**).

The autoclave validation is done in three steps: loggers programming, positioning inside the autoclave, reading of the data recorded during the sterilization cycle.

**Logger programming**

- Creation of a customized sterilization cycles library with spotting of the location of the loggers in use.

**Loggers programming**

- for an acquisition duration
- with an acquisition rate
- for a start at a specific date or on temperature threshold

## 1. Positioning of the loggers in autoclave

Placed inside the charges, TMI ORION probes can record any sterilization cycle type.

For average volume autoclaves, the basic model contains:


- 2 pressure/temperature loggers
- 4 temperature loggers
- A communication interface with the PC
- ValiVACQ data treatment software

For 1 cubic meter or more autoclaves, we recommend:

- 6 additional temperature embedded loggers
- 1 additional communication interface that connects to the first multi-tool interface; thus enabling simultaneous programming and reading of 12 loggers.

## 2. Simultaneous reading of recorded data

When the sterilization cycle is over, reconnect the loggers on the multi-tool interface. A few seconds later, all the data recorded in the NanoVACQ memory are displayed on your screen.

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Data analysis report in compliance with EN554 norm.

### System requirements

- Windows® 98-2000-Me-NT-XP,
- Serial port (RS232) or USB,
- Minimum 800x600 pixels display recommended